³He Top Loading Cryostat User Manual

Table of Contents:

Introduction:	. 2
Valve Names:	. 3
Sample Loading:	. 3
Liquid Nitrogen Cooldown:	. 4
Liquid Helium Cooldown	. 5
Cooling to Base Temperature:	
Sample Change:	. 8
. Temperature Contol	. 9
System Warm-up	10
Drawings & Pictures	
•	Valve Names: Sample Loading: Liquid Nitrogen Cooldown: Liquid Helium Cooldown Cooling to Base Temperature: Sample Change: Temperature Contol System Warm-up

I. Introduction:

- This is a complicated system that should only be run by an experienced user or member of the sample environment team. This manual is not meant to be a substitute for that, but instead should be a point of reference and a learning tool. For more information contact: Evan Fitzgerald (301-975-6657) or Dan Dender (301-975-6225).
- 2. What makes this system unique is that it is a *top loading* ³He system. Along with this feature come some inherent problems that had to be overcome to get the system to work. To reach a base temperature of 300mK the sample stick must make solid contact with the ³He pot and, for thermal shielding it, must also make contact with the 1K pot and the ⁴He reservoir. To accomplish this a force-generator was designed which pushes down on the sample stick and mates matching surfaces on the stick and the inside of the sample tube. The mating points are then connected by copper braids to the respective parts of the system. The system, however, still needs more cooling power to over come the large thermal mass of the stick and assembly, so there are actually two identical, independent ³He systems (insert A and B).

II. Valve Names:

- 1. V-1A/B = Charcoal cooling control
- 2. V-2A/B = 1K pot needle valve
- 3. V-3 = Vacuum space safety pressure relief
- 4. $V-4A/B = {}^{3}He$ dump fill
- 5. V-5A/B = 1K pots pumping and over pressure
- 6. V-6A/B = Charcoal cooling flow gauge
- 7. V-7 = Air lock double o-ring evacuation
- 8. V-8 = Air lock evacuation
- 9. V-9 = Sample well gate valve
- 10. V-10 = Inner vacuum chamber (IVC) evacuation
- 11. V-11 = IVC safety pressure relief
- 12. V-12 = Vacuum jacket evacuation
- 13. V-13 = Helium reservoir exhaust/pressure relief

Note: See "Figure 2" for a labeled picture of the top of the cryostat, "Figure 3" for a labeled picture of the sorb cooling valves and helium reservoir exhaust/pressure relief valve and "Figure 5" for a labeled picture of the sample stick.

III. Sample Loading:

1. With the system at room temperature, mount the sample stick on top of the

cryostat and bolt it into place. If the IVC is known to be under vacuum connect

a turbo pump to V-8 and evacuate the air lock. If the IVC is at atmospheric

pressure leave the airlock alone and continue. It is a good idea to occasionally

clean and grease the sample stick, but the amount of grease should be very

small.

2. Open V9 and slide the stick slowly into the sample well until it is stopped by the 0.5K thermal anchor stage. Screw the force generator down until it offers resistance to compress the stick slightly. Do NOT tighten the force generator all the way. See "Figure 5" for a picture of the sample stick mounted on the cryostat.

- 3. After checking that the hose between V7 and V8 is properly connected close both valves and evacuate the IVC through V10 with a turbo pump. When the vacuum is less than 5 $\times 10^{-4}$ Torr disconnect the turbo pump and connect a leak detector to V10.
- 4. Leak check the stick's hoses, valves, seals and connectors. Using a rough pump evacuate the 1K pot through V-5A/B, leaving V-2A/B closed, then backfill with He gas while monitoring the leak detector. Leave the 1K pot filled with He gas and the IVC evacuated. Disconnect the leak detector.

IV. Liquid Nitrogen Cooldown:

- Add a small amount of ⁴He gas into the IVC through V-10 by attaching a rubber hose, pinching it off about 1 foot down and then briefly opening then closing the valve.
- 2. Using a rough pump evacuate the helium reservoir through V-13, backfill with He gas and repeat three times, the last time leaving it filled with He gas. Start filling the helium and nitrogen reservoirs with LN₂.
- 3. Attach the rough pump to V-5A/B and pump on the 1K pots. Throughout the fill occasionally open V-2A/B one turn and then close it. This will prevent the needle valve from getting stuck. Stop this procedure when the 1K pot reaches less than 100K or when you stop filling the He reservoir with LN₂.
- 4. Close V-2A/B and pump out the 1K pots then back fill with He gas and seal off.
- 5. Let the LN_2 in the nitrogen reservoir fill all the way to the top, but stop filling the helium reservoir after approximately one hour or when the LN_2 level meter

reads above 70%. This assumes the two are filling at around the same rate with the goal being to leave a few inches of LN_2 in the helium reservoir.

- 6. Gently tighten the force generator if necessary as the stick will shrink as it cools, but do not over-tighten.
- 7. Let the system sit overnight to completely cool to Nitrogen temperature.

V. Liquid Helium Cooldown

- 1. If by the next morning the sensors read greater than 100K add more nitrogen to the helium reservoir and let things cool completely.
- 2. If the sensors are all reading near LN₂ temperature blow the LN₂ out of the helium reservoir by inserting a tube into the fill port and applying 3-5psi with helium gas into V-13. The fill port can slide up and down, so it needs to be pushed all the way down and then screwed into place. The nitrogen will flow out the tube, so use a rubber hose to direct it into a cryogenic liquid container.
- When all of the LN₂ is out use a rough pump to evacuate the helium reservoir through V-13 with V-2A/B open, backfill with He gas and repeat three times. Leave the reservoir filled with He gas the last time and then close V-2A/B.
- 4. Start transferring liquid helium (LHe) through the fill port after unscrewing it and slightly raising it. Remember to keep the pressure low (<3psi) and to remove V-13 to allow the exhaust to exit. See "Figure 3" for a picture of the helium exhaust valve.
- The LHe is full when the level meter reads >100% and a large dense plume of vapor exhausts through V-13. Remove the transfer line then seal the reservoir.
 V-13 should be sealed with the over-pressure valve to maintain a pressure of

about 3psi in the helium reservoir. See "Figure 3" for a picture of the helium exhaust valve.

- 6. Monitor the temperatures of both ³He inserts to make sure they cool steadily. When their temperatures reach approximately 25K start evacuating the IVC through V-10 using a turbo pump and allow at least 3 hours of pumping time. Do not allow the system to pump out longer than 8 hours as the sample temperature will slowly rise and be unacceptably high.
- 7. Without the exchange gas in the IVC cooling of the sample stick is dependent on the contact to the multiple mating surfaces. The force generator should be screwed tight, but not excessively so, to maintain contact with the cryostat.

VI. Cooling to Base Temperature:

- If the helium was recently filled wait a few minutes for some pressure to build up in the reservoir. Open V-1 and V-6 and verify flow by seeing the small black indicator ball jump up to a flow rate. If the ball doesn't move, even with a light tap, close V-1 immediately. The absence of flow is a problem that should be promptly reported to the sample environment team. If flow is verified close both V-1 and V-6 and continue on. See "Figure 3" for a picture of the valves involved in this procedure.
- If you are certain all of the He exchange gas is out of the IVC begin pumping on the 1K pot. Connect a high volume rough pump to V-5A/B and pump out the lines back to V-2A/B. Slowly open V-2A/B and adjust to reach the lowest sustainable base temperature as measured on the 1K pot sensors (~1.1 – 1.3 K).
 See "Figure 4" for a picture and description of the electronics rack setup.

- 3. When the 1K pots are steadily at their base temperatures set the temperature of the ³He sorbs directly to 20K (see below to turn the ramp feature off) and turn their heater on to "Medium". See "Figure 4" for a picture and description of the electronics rack setup. The 1K pot temperatures will briefly rise while they cool the "warm" ³He gas. After the sorb has reached 20K turn on the ramp feature by pressing the "control setup" button on the sorb controller and using the "enter" button to scroll to the ramp option. Set the ramp on and then set the rate to 1 K/m. With the ramp enabled change the setpoint to 40K and the controller should slowly ramp there. The sample temperature should slowly decrease to around 1.8K. The ³He is condensing when the ³He pots reach 1.8K. Continue to condense for at least half an hour after this point (1 hour condensation = ~25hours at base; 4 hours of condensation = ~39hours at base).
- 4. Once the ³He has condensed for the desired amount of time turn off the charcoal adsorption pump heater and slowly open V-1 and V-6. Once again watch to see the black ball on the flow meter measure a flow rate and close the valve if it doesn't move. The recommended flow rate is approximately 1,000cc per minute and the charcoal adsorption pump should slowly cool to, and remain below, 5K.
- 5. The ³He pot temperatures will slowly drop to 300mK followed closely by the sample temperature. The 1K pots are the main contributors to helium consumption when running and should be monitored and adjusted as necessary to lower their temperature and minimize this consumption. The Silicon Diode type sample sensor (A), which bottoms out at 1.4K should be disabled to remove

7

the heat load it contributes. Press the "Input Setup" button on the sample sensor's controller and change the first setting, "Enable", from "on " to "off". See "Figure 4" for a picture and description of the electronics rack setup.

VII. Sample Change:

- Connect a turbo pump to V-8, leaving it closed, and evacuate the air lock assembly. See "Figure 5" for a picture of the airlock and associated valves.
 Open V-7 and begin to slowly raise the sample stick. As the stick is pulled up it will frost up and can freeze the o-ring seal on the air lock. To avoid this pause every few minutes and allow the stick to warm or heat it with a heat gun. The clamp can be used to hold the stick in place as you are warming it up.
- 2. When the stick is all the way up and will not go any farther clamp it into place, close V-9 and V-7 and stop the pump. The air lock and stick can be brought to atmospheric pressure by slowly opening V-8. Helium gas can alternatively be added to the air lock to slowly warm the sample if necessary. (Note: Helium gas will expand as it warms and will quickly build up pressure in a trapped volume like this. Contact a member of the sample environment team before attempting this procedure.)
- 3. Unbolt the stick from the top of the cryostat, remove it and cover V-9 with the blank flange to protect it.
- To load another sample re-mount the stick and bolt it down. Evacuate the air lock with V-7 and V-8 open.
- Have an experienced user add a small amount of ³He gas to the IVC through V-10, close V-8, open V-9 and begin slowly lowering the stick. The use of ³He

8

instead of ⁴He will decrease the amount of time required to pump out the IVC and should only be used when time is a major factor. Continue pumping on V-7 to prevent any air from getting by the double o-ring seal. The LHe will begin to rapidly boil off due to the large heat load. To decrease the helium boil off lock the stick in place and let it sit for fifteen minutes every few feet. The stick should be as cold as possible, but certainly less then 100K, before making contact with the mating surfaces in the sample well.

- When the stick sensors reach approximately 25K pump out the IVC through V-10 using a turbo pump for at least a few hours.
- 7. Continue following the directions from step 5 of section V.

VIII. Temperature Contol

- T < 1K: When starting from base temperature the sorb can be gently warmed (T_{sorb} < 15K, reduce flow through V6) to increase sample temperature slightly. This method can actually increase the time you can stay below 1.8K. Using the sample heater at its lowest power setting is another option, but this will dramatically reduce the amount of time you can stay below 1.8K.
- 2. 1 K < T > 1.8 K: The ³He must be condensed and pumped on to work below 1.8K, but it is difficult to make it last for more than a few hours when working above 1K. To work in this temperature range use both the sample heater and the sorb heater. Warming the sorb to about 20 or 30K will slow the pumping of the ³He and then the sample heater can be used at a low power to control the temperature of the sample.

- 3. 1.8K < T > 4.2K: In this range cooling from the 1K pots is necessary, but pumping on liquid ³He is not. The sorb temperature should be set at 25K or above to provide exchange gas in the ³He circuit. The sample heater can then be used at an appropriate power.
- 4. 4.2K < T > 300K: At temperatures close to 4 K the same method as suggested in #3 should be employed. At higher temperatures the cooling from the 1K pot will become useless and they can be closed off to conserve LHe. At this point the sorb temperature is also not important, so its heater can be turned off. The sample will then be isolated and its temperature entirely dependant on the heater and radiative cooling. Above 20K the RuOx sample sensor (B) will not read, so the control sensor in the 340 controller and the sample temperature sensor read by ICP must be changed to the Silicon Diode sensor (A).
- Note: For temperatures above approximately 30K it is necessary to add exchange gas to cool the sample back down at a reasonable rate. This exchange gas will then take a few hours to pump back out if you wish to continue below 4K.

IX. System Warm-up

- Valves 2A/B and 5A/B should be closed to stop pumping on the 1K pots. Although there is cold He trapped between valve 2 and 5 there is an overpressure valve that will safely release any excess pressure.
- 2. Close valves 1A/B to seal the cooling circuit for the charcoal pumps and prevent air from being sucked in when the LHe runs out.

- 3. Finally loosen the force generator slightly so that it will be easier to remove at room temperature after the stick has expanded.
- 4. As long as the LHe reservoir is sealed as normal the system should be safe to warm up on its own.

X. Drawings & Pictures




Figure 1: Cryostat Side View Drawing with Dimensions.




Figure 2: Cryostat Top View With Parts Labeled




Figure 3: Sorb Cooling Circuit Operation and Filling Liquid Helium






Figure 5: Sample Stick and Air Lock Parts