

## Top Loading CCR PID Settings

For low temperature sample stick with 1psi of helium in sample well.

Zone #	Zone Limit (K)	P	I	D	Manual Output	Power
1	55	150	10	0	0	50W
2	150	150	10	20	0	50W
3	200	100	20	30	0	50W
4	240	75	20	30	0	50W
5	280	60	20	30	0	50W

For high temperature sample stick with vacuum in sample well.

Zone #	Zone Limit (K)	P	I	D	Manual Output	Power
1	375	200	2	8	0	50W
2	425	240	0.1	8	0	50W
3	725	240	0.2	8	0	50W

9/29/2006